

Tested cleaning products

- Always use the disinfectant according to the instructions of the manufacture.
- We recommend to test a disinfectant on an invisible point of the covering before any regular application of this disinfectants.

MANUFACTURER	DISINFECTION
Akzo Nobel	Natriumhypochlorid
Bode Chemie GmbH	Bacillol AF Bacillol 30 Sensitive Foam Bacillol 30 Sensitive Tissues Cutasept Mikrobac forte
Braun	Meliseptol Foam Pure Softasept N Meliseptol Helipur Melsitt Melsept
Desomed	Desomed Rapid
Dr. Schumacher	Descosept AF
Ecolab	Incidin foam Incidin Plus Incidin perfekt Incidin liquid
Favodent Karl Huber	Flora quick
Henkel	Incidur-Spray
Lysoform	Fugaten Spray
Orochemie	B45
Peppler	Micro clean Pepp clean
Ruck	Sprühdesinfektion
Reckitt Benckiser	Sagrotan med
Schülke & Mayr	Kodan Tinktur Forte Esemfix Octenisept Mikrozid AF
Unident	Unisepta foam